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Issuing organization FOI – Swedish Defence Research Agency NBC Defence SE-901 82 Umeå	Report number, ISRN FOI-R--1637--SE	Report type Methodology report
	Programme Areas 3. NBC Defence and other hazardous substances	
	Month year December 2005	Project no. A 438
	Subcategories 32 Biological and Chemical Defence Research	
	Subcategories 2	
Author/s (editor/s) Daniel Johansson Victoria Heldestad Ann Göransson Nyberg	Project manager Ann Göransson Nyberg	
	Approved by	
	Sponsoring agency Department of defence	
	Scientifically and technically responsible	
Report title (In translation) Method for blood flow measurements with fluorescent microspheres		
Abstract <p>The principle of blood flow measurements using microspheres is simple and ingenious. By injecting microspheres directly into the heart of experimental animals, the spheres will be spread through blood circulation. Since the spheres have a larger diameter than the smallest capillaries in the different tissues of the body, they will be captured when trying to get through the capillary bed. Assuming that microspheres are mixed homogenously when they leave the heart, the spheres that are trapped within different tissues will be in proportional relation to the blood flow in that tissue at the time of the injection. Previously, radioactive spheres have been used to measure how much is trapped in the tissue. It is possible to measure blood flow at different times in the same animal, if different isotopes are used. The problems associated with radioactive microspheres are short half-life, high costs and risks for the users. Today, different fluorescent microspheres can be used as markers with advantages such as long validity, low costs and safety. The disadvantage is that the spheres have to be prepared from the tissue prior to quantification. The purpose of the present study is to establish a method to measure blood flow using fluorescent microspheres, in experimental animals, at FOI NBC Defence in Umeå.</p>		
Keywords Fluorescent microspheres, blood flow, cerebral, regional		
Further bibliographic information	Language Swedish	
ISSN 1650-1942	Pages 27 p.	
	Price acc. to pricelist	

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Referenser

de Garavilla L, Dalsey W, Clas D, Eynon CA. Inter- and intra-animal variability in baseline organ blood flow in rats as measured with fluorescent microspheres.

Abstract, Third International Conference on fluorescent microsphere methods 1996, Seattle, Washington.

el-Etri MM, Nickell WT, Ennis M, Skau KA, Shipley MT. (1992). Brain norepinephrine reductions in soman-intoxicated rats: association with convulsions and AChE inhibition, time course, and relation to other monoamines. *Exp Neurol.* Nov;**118**(2):153-63.

Elhusseiny A, Cohen Z, Olivier A, Stanimirovic' DB, Hamel E. (1999). Functional Acetylcholine Muscarinic Receptor Subtypes in Human Brain Microcirculation: Identification and Cellular Localization. *J Cereb Blood Flow Metab* **19**: 794-802.

Glenny, R. W., S. Bernard and M. Brinkley.(1993). Validation of fluorescent-labeled microspheres for measurement of regional organ perfusion. *J Appl Physiol.* **74**:2585-97.

Göransson Nyberg A, and Cassel G. (2001). Cardiopulmonary effects of HI 6 treatment in Soman intoxication. *J Appl Toxicol*, **21**, 79-81.

Koskinen, L-O.D. (1989). Cerebral and peripheral blood flow effects of TRH in the rat-A role of vagal nerves. *Peptides.* **10**, 933-938.

Manual for Using Fluorescent Microspheres to Measure Regional Organ Perfusion

<http://fmrc.pulmcc.washington.edu/fmrc.shtml>

Maxwell, D.M., Lenz, D.E., Groff, W.A., Kaminskis, A., Froehlich, H.L., (1987). The effects of blood flow and detoxification on in vivo cholinesterase inhibition by soman in rats. *Toxicol. Appl. Pharmacol.* **88**, 66-76.

McDonough JH, Shih TM. (1993). Pharmacological modulation of soman-induced seizures. *Neurosci Biobehav Rev* **17**:203-215.

Prinzen, F. W. and R. W. Glenny. (1994). Developments in non-radioactive microsphere techniques for blood flow measurement. *Cardiovasc Res.* **28**:1467-75.

Takeuchi Y, Matsushita H, Kawano H, Sakai H, Yoshimoto K, Sawada T. (1999). TRH increases cerebrospinal fluid concentration of kynurenine. *Neuroreport.* Nov **26**, 10(17):3601-3.

Taylor, P. Anticholinesterase agents. In: The Pharmacological basis of therapeutics Goodman Gilman ed. by A, Rall TW, Nied AS, Taylor P Pergamon Press, New York (1990).

van Oosterhout, M. F., H. M. Willigers, R. S. Reneman and F. W. Prinzen.(1993). Fluorescent microspheres to measure organ perfusion: validation of a simplified sample processing technique. *Am J Physiol.* **269**:H725-33.

Wagner HN Jr, Rhodes BA, Sasaki Y, et al: (1969). Studies of the circulation with radioactive microspheres. *Invest Radiol* **4**:374-386.

Bilaga I: Description of the Sample Preparation tube set

Introduction

The Sample Processing tube system, which was designed and validated by the Institute for Surgical Research at the Uni-Klinik Großhadern in Munich, contains a set of parts (Figure 1) which can be used for all of the sample preparation steps for blood flow measurement.


Figure 1. The Sample Preparation tube set (B0511133) consists of the filter (A), the filter holder (B), screw-top cap for the filter holder/centrifuge tube (C), screw-top cap for the measurement tube (D) and measurement tube (E).

Figure 2. The stainless steel digestion tube (F) and delrin cap (G).

Also delivered with the Blood flow system (but not included in the Sample Preparation tube set) are the stainless steel sample preparation tube (F B801-8081, included with the heater block) and the delrin cover (G, B801-8098, pack of 40):

The central part of the set is the filter, which has a polyamide membrane permanently fixed into its base.

This membrane allows the vacuum filtration of digested tissue samples and arterial reference samples. The 7 micron pore size of the filter guarantees the complete retention of all of the 15 micron fluorescence microspheres which are extracted from the tissue and blood samples.

Tissue samples remain in the same filter tube throughout the complete sample preparation procedure, eliminating the possibility of loss of microspheres.

Preparation protocol for tissue samples

Reagents

Aqueous potassium hydroxide (4N)

222,44 g KOH pellets are dissolved in 1000 ml deionised water, this solution is stirred and 10 ml 2% Tween 80 added.

Isopropanol (100% !!!)

Phosphate buffer

29.90 g dipotassium hydrogen phosphate trihydrate is dissolved in 800 ml deionised water. 5,88 g potassium hydrogen phosphate is dissolved in 200 ml deionised water.

The two solutions are mixed together.

Method

After weighing, each tissue sample is added to a filter tube. The filter tube is then inserted into the stainless steel sample preparation tube:


Figure 3. Insertion of the filter (A) into the stainless steel digestion tube (F).

Figure 4. The delrin cover (G) placed onto the stainless steel digestion tube (F).

15 ml of aqueous potassium hydroxide reagent is then pipetted into the filter tube. To this mixture 1.5 ml isopropanol (100% !!!) is added, then the delrin cover is placed over the tube to minimise evaporation during digestion:

Digestion is carried out in the heater block (B801-8064) at 60C. Digestion of tissue samples takes 6 hours to completion, although samples with a weight of less than 3 g are typically completely dissolved within 4 hours.

The chemical and physical properties of the Sample Preparation tube set also allow the preparation of strongly lipophilic tissue such as brain.

Following digestion, the filter is removed from the stainless steel digestion tube and placed onto the vacuum station (B801-8061). After a few seconds, the KOH is removed by the action of the vacuum. All remaining traces of KOH are then removed by thoroughly washing the filter with 30 ml phosphate buffer, which is applied in a circular motion over the whole surface and round the edges of the filter. This is done with the filter's in and outside.

The filter is then centrifuged at 2.800 r.p.m. for 1-2 minutes to ensure that the filter is dry.

Note that if the filter is not thoroughly dry, then problems may occur with solution of the spheres and collection of the fluorescent marker solution, due to the immiscibility of water and the solvent. The filter tube can be placed directly in a polypropylene centrifuge tube. Using the filter holder cap will prevent leakage of (corrosive) KOH into the centrifuge:


Figure 5. Inserting the filter (A) into a polypropylene centrifuge tube (H) and fixing the screw cap (C). Note that the standard Perkin-Elmer 50 ml polypropylene tube and screw-top cap (B0193234) are compatible with the filter system.

Figure 6. The filter (A) is inserted into the holder (B) and the measurement tube (D) attached. 2 ml of solvent is added to the filter tube, solution of the microspheres is then assisted by mechanical shaking and vortexing gently for 1-2 minutes. The solvent is collected in the measurement tube by centrifuging the whole tube system at 2,800 r. p. m. for 2 minutes. It is then ready for placing onto the auto sampler tray for measurement.

Figure 7. Note that if it is not possible to carry out measurement immediately, then the cap must be fixed onto the tube.

Failure to do this could lead to evaporation of solvent, which causes an increase in dye concentration and correspondingly false results.

The filter is then removed from the centrifuge tube and inserted into the filter holder, the measurement tube is then fixed to the filter holder.

Preparation protocol for lung tissue samples

Reagents

Ethanolic potassium hydroxide (2N).

112.22 g KOH pellets are dissolved in 1000 ml 96% ethanol.

Phosphate buffer (see above).

(No isopropanol !!!)

Methods

The preparation protocol for lung samples is similar to that described above, except that ethanolic KOH is used instead of aqueous KOH, and the addition of isopropanol must not take place.

Preparation protocol for bone tissue samples

Reagents

1N HCl.

Aqueous potassium hydroxide (4N, see above).

Isopropanol (100% !!!), 3 ml

Phosphate buffer (see above).

Methods

Preparation of bone samples differs from other samples in that the crystalline matrix of the bone must be dissolved with 1N HCl before digestion with KOH is carried out.

Bone samples are inserted into a polypropylene tube and are treated with acid for 1 week followed by washing with phosphate buffer. The samples are then inserted into the filters. Thereafter the method is identical to the procedure outlined in “(2) Preparation protocol for tissue samples”.

Preparation protocol for blood samples

Reagents

Isotonic saline, 0.9% w/v, containing 0.5% v/v tween 80

CPDA-1 anticoagulation stabilisor.

26,3 g sodium citrate dehydrate, 3.27 g citric acid monohydrate, 31.9 g glucose monohydrate, 2.51 g sodium dehydrogenate phosphate dehydrate and 0.275 g adenine are dissolved in 1000 ml water for injection.

OR!!! Na-Citrate: 1 ml of Na-Citrate for 2 ml of blood

Alternatively the reference blood samples can be anticoagulated by adding Citrate (1 part Citrate : 2 parts blood) to the collecting syringe prior to blood withdrawal.

Aqueous potassium hydroxide (4N, see above).

Isopropanol (100% !!!)

Phosphate buffer (see above).

Methods

The syringe to be used for withdrawal of the reference blood sample must contain CPDA-1 anticoagulation stabilisor, 1:5 v/v stabilisor: blood, prior to withdrawal of the reference sample.

After withdrawal of the sample, the contents of the entire withdrawal system (syringe and catheter) are filtered through the filter system using the vacuum station. Following this, the withdrawal system is washed with excess isotonic saline or phosphate buffer (see above), and the washing solution is filtered again through the filter using the vacuum station.

The samples are then treated as described above in “(2) Preparation protocol for tissue samples”.

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Construction of Standard Curves

- 1) Take 5 μl of 0.2% FluoSpheres® (5,000 fluorescent microspheres) and add it to 10 mls of solvent (we use 2-(2-ethoxyethoxy) ethyl acetate). This yields the fluorescence intensities for ~ 500 fluorescent microspheres/ml. Make an emission scan with a fixed excitation wavelength (use the wavelength recommended on the bottle). Note the best emission wavelength and perform an excitation scan with a fixed optimal emission wavelength. Make sure that the excitation peak and emission peak can be distinguished and do not interfere with the other colors used. Play around with the bandpasses to get a decent strength of your signal.
- 2) Depending on the range of your fluorimeter (based on the manufacturer's recommended operating range), calculate the number of μl required to add to 10 mls solvent to yield a fluorescent microsphere concentration that will be just above the upper intensity limit for your fluorimeter.
- 3) Using this high-stock solution, make multiple serial dilutions of microspheres with very accurate pipetting. Calculate the number of microspheres per ml for each dilution.
- 4) The fluorescence intensities per number of microspheres per ml are plotted to yield a standard curve for that lot of microspheres.

Making Control Solutions

When making an appropriate set of control solutions, it is important to consider the number of fluorescent microspheres per volume of solvent being used and their "typical" intensities.

- 1) Control solutions should have a fluorescent signal that is approximately the mean signal for each color in your sample, and be on the linear part of the standard curve, as described earlier.
- 2) Take 5 μl of 0.2% FluoSpheres® (5,000 fluorescent microspheres) and add it to 10 mls of solvent (we use 2-ethoxyethyl acetate). This yields fluorescence intensities of ~ 500 fluorescent microspheres/ml and gives a reference intensity for making your control solution.

3) Depending on the number of different colored microspheres being used ("n"), a stock solution is made by dissolving one color of microspheres in 2-ethoxyethyl acetate to yield an intensity that is the "n" times more concentrated than average intensities from an experimental sample.

4) Stock solutions are then combined to yield a control solution in the range of our "typical" sample intensities.

For example, if we use three fluorescent microsphere colors (blue-green, orange and crimson), and our "typical intensities" are 100 for blue-green, 200 for orange and 50 for crimson, we should make a stock solution 3 x 100, or 300 for blue-green, 3 x 200, or 600 for orange, and 3 x 50, or 150 for crimson. The three stock solutions are then combined to yield a 1:3 dilution to make a control solution that reads 100, 200 and 50, respectively, for blue-green, orange and crimson.

Method for Reference Blood Flow Sampling

There are 2 different anticoagulants that we routinely use: Heparin (syringe coated) and Citrate Phosphate Dextrose (10 cc per 30 cc syringe).

1. Using whole blood, calibrate the reference withdrawal pump at the predetermined withdrawal rate, including the catheters, extension tubing and matched anticoagulated glass syringes that will be used for the reference withdrawal. Have new stopcocks and flush syringes available.
2. Connect the matched glass anticoagulated syringes in the withdrawal pump to the catheters and the extension tubing so that everything is set up for withdrawing the reference sample. Do not turn the stopcock on the catheters until you are ready for injection (the blood will flux into the catheter dead space and may clot).
3. Once the microspheres have been drawn into the injection syringes, start the withdrawal pump and make sure blood is flowing freely into the extension tubing.
4. Inject the microspheres over the designated time period (sec or min) followed by a flush of warmed saline three times the volume of the catheter dead space.
5. A timer is set for 2 min after completion of injection for the reference blood withdrawal. At the end of the withdrawal, the pump is turned off, the stopcocks are opened and the blood remaining in the extension tubing is drawn into the syringe.
6. Transfer blood into labelled vials for further processing (see Digestion of Blood and Tissue, next page). Rinse syringes and extension lines with 2% Tween-80® (using approximately twice the volume of the blood) and add this rinse to the blood samples.
7. Flush the catheters again and change the stopcocks.